



KIR2DL5A Rabbit pAb

Catalog No	YP-Ab-17731
Isotype	IgG
Reactivity	Human
Applications	WB,IHC-P
Gene Name	KIR2DL5A
Alternative Names	killer cell immunoglobulin like receptor; two Ig domains and long cytoplasmic tail 5A; CD158F; KIR2DL5; KIR2DL5.1; KIR2DL5.3
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa
Clonality	Polyclonal Antibody
Dilution	WB: 1/500-1/1000 IHC: 1/50-1/100
Immunogen	Synthetic peptide of human KIR2DL5A
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Receptor on natural killer (NK) cells for HLA-C alleles.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

