



ApoA2 Rabbit mAb

Catalog No	YP-Ab-17769
Isotype	IgG
Reactivity	Human,Rat
Applications	WB,ICC/IF,FC
Gene Name	APOA2
Alternative Names	APO A2; Apo AII; APOA 2; ApoA II; APOA2; APOAII; Apolipoprotein A II
Research Field	Cardiovascular
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa
Clonality	Monoclonal Antibody
Clonality No.	R09-8W8
Dilution	WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100
Immunogen	A synthesized peptide derived from human Apolipoprotein A II
Purification	Affinity Chromatography
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

