



## MTX2 rabbit pAb

货号	YP-Ab-17793
同位型	IgG
种属	Human;Mouse
应用	WB
靶点	MTX2
基因名称	MTX2
蛋白名称	Metaxin-2 (Mitochondrial outer membrane import complex protein 2)
免疫原	Synthesized peptide derived from human MTX2
特异性	This antibody detects endogenous levels of MTX2 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	29kD
功能	Involved in transport of proteins into the mitochondrion.
细胞定位	Mitochondrion outer membrane . Mitochondrion .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images