



TRIM29 Rabbit mAb

Catalog No	YP-Ab-17819
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TRIM29
Alternative Names	ATDC; TRIM29
Research Field	Epigenetics and Nuclear Signaling
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 66 kDa; Observed MW: 66 kDa
Clonality	Monoclonal Antibody
Clonality No.	R02-3S5
Dilution	WB: 1/500-1/1000
Immunogen	A synthesized peptide derived from human TRIM29
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	It is able to complement the radiosensitivity defect of an ataxia telangiectasia (AT) fibroblast cell line.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

