

## Mouse Anti Human sFLT-1 mAb-Clone 2

**Catalog No:**

**Quantity:** Specified on product label

**Buffer:** PBS,PH7.4

**Storage:** Store at 2-8C short term and -80°C long term.Storage

It is recommended to aliquot needed volumes and store at -80C  
Avoid repeated freeze/thaw cycles

### QC Results

| Characteristics | Method                                   | Requirements | Results   |
|-----------------|--|--------------|-----------|
| Appearance      | Visual observation                       | Clear Liquid | Qualified |
| Concentration   | A280 with the extinctioncoefficient 1.43 | >2.00mg/mL   | 8.26mg/mL |
| Purity          | SDS-PAGE                                 | 295%         | >95%      |