

FOXO6 rabbit pAb

货号	YP-Ab-17835
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	FOXO6
基因名称	FOXO6
蛋白名称	Forkhead box protein O6
免疫原	Synthesized peptide derived from human FOXO6
特异性	This antibody detects endogenous levels of FOXO6 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	54kD
功能	Transcriptional activator.
细胞定位	Cytoplasm . Nucleus . When phosphorylated, translocated from nucleus to cytoplasm. High nuclear localization after stimulation with growth factors (By similarity). .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images