

RPS29 rabbit pAb

货号	YP-Ab-17836
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	RPS29
基因名称	RPS29
蛋白名称	40S ribosomal protein S29
免疫原	Synthesized peptide derived from human RPS29
特异性	This antibody detects endogenous levels of RPS29 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	6kD
细胞定位	Cytoplasm, cytosol . Cytoplasm . Rough endoplasmic reticulum . Detected on cytosolic polysomes (PubMed:25957688). Detected in ribosomes that are associated with the rough endoplasmic reticulum (By similarity). .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images