



# CCR6 Rabbit pAb

<b>Catalog No</b>	YP-Ab-17863
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human
<b>Applications</b>	WB,IHC-P
<b>Gene Name</b>	CCR6
<b>Alternative Names</b>	CCR6; CKRL3; CMKBR6; GPR29; STRL22; C-C chemokine receptor type 6; C-C CKR-6; CC-CKR-6; CCR-6; Chemokine receptor-like 3; CKR-L3; DRY6; G-protein coupled receptor 29; GPR-CY4; GPRCY4; LARC receptor; CD196
<b>Research Field</b>	Immunology
<b>Product Categories</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Molecular Weight</b>	Calculated MW: 42 kDa; Observed MW: 42 kDa
<b>Clonality</b>	Polyclonal Antibody
<b>Clonality No.</b>	-
<b>Dilution</b>	WB: 1/500-1/1000 IHC: 1/50-1/100
<b>Immunogen</b>	Synthetic peptide of human CCR6
<b>Purification</b>	Affinity Purified
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Form</b>	Liquid
<b>Buffer System</b>	pH 7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Background</b>	Receptor for the C-C type chemokine CCL20.
<b>matters needing attention</b>	Avoid repeated freezing and thawing!



## Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images

