



# INPP4B Rabbit mAb

<b>Catalog No</b>	YP-Ab-17864
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human,Rat
<b>Applications</b>	WB
<b>Gene Name</b>	INPP4B
<b>Alternative Names</b>	4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II 105kDa; Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3; 4 bisphosphate 4 phosphatase; Type II inositol-3
<b>Research Field</b>	Signal Transduction
<b>Product Categories</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Molecular Weight</b>	Calculated MW: 105 kDa; Observed MW: 105 kDa
<b>Clonality</b>	Monoclonal Antibody
<b>Clonality No.</b>	R08-7C7
<b>Dilution</b>	WB: 1/500-1/1000
<b>Immunogen</b>	A synthetic peptide of human INPP4B
<b>Purification</b>	Affinity Purified
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodifi
<b>Form</b>	Liquid
<b>Buffer System</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Background</b>	4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II 105kDa; Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3; 4 bisphosphate 4 phosphatase; Type II inositol-3
<b>matters needing attention</b>	Avoid repeated freezing and thawing!



## Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images

