



FNDC5 rabbit pAb

货号	YP-Ab-17883
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	FNDC5
基因名称	FNDC5 FRCP2
蛋白名称	Fibronectin type III domain-containing protein 5 (Fibronectin type III repeat-containing protein 2) [Cleaved into: Irisin]
免疫原	Synthesized peptide derived from human FNDC5
特异性	This antibody detects endogenous levels of FNDC5 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	23kD
功能	[Irisin]: Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of browning of human white adipose tissue.
细胞定位	Cell membrane ; Single-pass type I membrane protein . Peroxisome membrane ; Single-pass type I membrane protein . Imported in peroxisomes through the PEX5 receptor pathway. .; [Irisin]: Secreted. Detected in the blood of individuals subjected to endurance exercise. .
组织表达	Widely expressed, with highest levels in heart. Very low expression, if any, in colon, pancreas and spleen.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

