



# HOXD10 Rabbit pAb

<b>Catalog No</b>	YP-Ab-17886
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human
<b>Applications</b>	WB,IHC-P
<b>Gene Name</b>	HOXD10
<b>Alternative Names</b>	HOX4; HOX4D; HOX4E; Hox-4.4
<b>Research Field</b>	Epigenetics and Nuclear Signaling
<b>Product Categories</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Molecular Weight</b>	Calculated MW: 38 kDa; Observed MW: 38 kDa
<b>Clonality</b>	Polyclonal Antibody
<b>Clonality No.</b>	-
<b>Dilution</b>	WB: 1/500-1/1000 IHC: 1/50-1/100
<b>Immunogen</b>	Synthetic peptide of human HOX80
<b>Purification</b>	Affinity Purified
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Form</b>	Liquid
<b>Buffer System</b>	pH 7.4 PBS, 0.05% NaN3, 40% Glycerol
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Background</b>	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
<b>matters needing attention</b>	Avoid repeated freezing and thawing!



### Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images

