



DEPDC1 Rabbit pAb

Catalog No	YP-Ab-17905
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	DEPDC1
Alternative Names	DEP domain containing 1; DEP.8; SDP35; DEPDC1A; DEPDC1-V2
Research Field	-
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 93 kDa; Observed MW: 93 kDa
Clonality	Polyclonal Antibody
Clonality No.	-
Dilution	WB: 1/500-1/1000
Immunogen	Fusion protein of human DEPDC1
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	May be involved in transcriptional regulation as a transcriptional corepressor.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

