



Hippocalcin Rabbit mAb

Catalog No	YP-Ab-17907
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HPCA
Alternative Names	BDR2; HPCA; P23K
Research Field	Neuroscience
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 22 kDa; Observed MW: 26 kDa
Clonality	Monoclonal Antibody
Clonality No.	R08-1N7
Dilution	WB: 1/500-1/1000
Immunogen	A synthesized peptide derived from human Hippocalcin
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	May be involved in the calcium-dependent regulation of rhodopsin phosphorylation. Binds two calcium ions.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

