



Phospho-RIP3 (Ser232) Rabbit mAb

Catalog No	YP-Ab-17922
Isotype	IgG
Reactivity	Mouse
Applications	WB
Gene Name	RIPK3/Ripk3
Alternative Names	Receptor interacting protein 3; Receptor interacting serine threonine kinase 3; Receptor interacting serine/threonine protein kinase 3; Receptor-interacting serine/threonine-protein kinase 3; RIP like protein kinase 3; RIP-3;RIP-like protein kinase 3; RIPK3;
Research Field	Cell Biology
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW:53/57kDa;Observed MW:50kDa
Clonality	Monoclonal Antibody
Clonality No.	R66-4K-6
Dilution	WB: 1/500-1/1000
Immunogen	A synthesized peptide derived from human Phospho-RIP3 (S232)
Purification	Affinity Chromatography
Conjugation	Unconjugated
Modification	Phosphorylated
Form	Liquid
Buffer System	Liquid in PBS , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Promotes apoptosis.
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

