

AXUD1 rabbit pAb

货号	YP-Ab-17944
同位型	IgG
应用	WB
种属	HHuman;Mouse;Rat
靶点	MIOX
基因名称	MIOX ALDRL6 KSP32 RSOR
蛋白名称	Inositol oxygenase (EC 1.13.99.1) (Aldehyde reductase-like 6) (Kidney-specific protein 32) (Myo-inositol oxygenase) (MI oxygenase) (Renal-specific oxidoreductase)
免疫原	Synthesized peptide derived from human MIOX
特异性	This antibody detects endogenous levels of MIOX at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	31kD
细胞定位	Cytoplasm .
组织表达	Kidney specific.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images