

RAB3D rabbit pAb

货号	YP-Ab-17952
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	RAB3D
基因名称	RAB3D GOV RAB16
蛋白名称	Ras-related protein Rab-3D
免疫原	Synthesized peptide derived from human RAB3D
特异性	This antibody detects endogenous levels of RAB3D at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	24kD
功能	Protein transport. Probably involved in regulated exocytosis (By similarity).
细胞定位	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
组织表达	Highly expressed in granulocytes of peripheral blood. Constitutively expressed at low levels in all hematopoietic cell lines investigated.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images