



AXIN2 rabbit pAb

货号	YP-Ab-17954
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	AXIN2
基因名称	AXIN2
蛋白名称	Axin-2 (Axin-like protein) (Axil) (Axis inhibition protein 2) (Conductin)
免疫原	Synthesized peptide derived from human AXIN2
特异性	This antibody detects endogenous levels of AXIN2 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	93kD
功能	Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B.
细胞定位	Cytoplasm .
组织表达	Expressed in brain and lymphoblast.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images