

KLHL14 rabbit pAb

货号	YP-Ab-17967
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	KLHL14
基因名称	KLHL14 KIAA1384
蛋白名称	Kelch-like protein 14 (Protein interactor of Torsin-1A) (Printor) (Protein interactor of torsinA)
免疫原	Synthesized peptide derived from human KLHL14
特异性	This antibody detects endogenous levels of KLHL14 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	69kD
细胞定位	Cytoplasm, cytosol . Endoplasmic reticulum membrane . Colocalizes with TOR1A at the endoplasmic reticulum level.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images