

TMX3 rabbit pAb

货号	YP-Ab-17975
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	TMX3
基因名称	TMX3 KIAA1830 TXNDC10
蛋白名称	Protein disulfide-isomerase TMX3 (EC 5.3.4.1) (Thioredoxin domain-containing protein 10) (Thioredoxin-related transmembrane protein 3)
免疫原	Synthesized peptide derived from human TMX3
特异性	This antibody detects endogenous levels of TMX3 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	50kD
功能	Probable disulfide isomerase, which participates in the folding of proteins containing disulfide bonds. May act as a dithiol oxidase.
细胞定位	Endoplasmic reticulum membrane ; Single-pass membrane protein .
组织表达	Widely expressed. Expressed in brain, testis, lung, skin, kidney, uterus, bone, stomach, liver, prostate, placenta, eye and muscle.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images