



HIGD1A rabbit pAb

货号	YP-Ab-17976
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	HIGD1A
基因名称	HIGD1A HIG1 HSPC010
蛋白名称	HIG1 domain family member 1A (Hypoxia-inducible gene 1 protein)
免疫原	Synthesized peptide derived from human HIGD1A
特异性	This antibody detects endogenous levels of HIGD1A at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	10kD
功能	Proposed subunit of cytochrome c oxidase (COX, complex IV), which is the terminal component of the mitochondrial respiratory chain that catalyzes the reduction of oxygen to water. May play a role in the assembly of respiratory supercomplexes.
细胞定位	Mitochondrion membrane ; Multi-pass membrane protein . Mitochondrion inner membrane .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images