



S3-12 rabbit pAb

货号	YP-Ab-17979
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	PLIN4
基因名称	PLIN4 KIAA1881
蛋白名称	Perilipin-4 (Adipocyte protein S3-12)
免疫原	Synthesized peptide derived from human S3-12
特异性	This antibody detects endogenous levels of S3-12 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	151kD
功能	May play a role in triacylglycerol packaging into adipocytes. May function as a coat protein involved in the biogenesis of lipid droplets (By similarity).
细胞定位	Cell membrane . Cytoplasm . Lipid droplet . Nascent lipid droplet surface-associated; association with lipid droplets is triacylglycerol synthesis-dependent. .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images