



UBE2C rabbit pAb

货号	YP-Ab-17997
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	UBE2C
基因名称	UBE2C UBCH10
蛋白名称	Ubiquitin-conjugating enzyme E2 C (EC 6.3.2.19) (Ubch10) (Ubiquitin carrier protein C) (Ubiquitin-protein ligase C)
免疫原	Synthesized peptide derived from human UBE2C
特异性	This antibody detects endogenous levels of UBE2C at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	20kD
功能	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-11'- and 'Lys-48'-linked polyubiquitination. Acts as an essential factor of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated ubiquitin ligase that controls progression through mitosis. Acts by initiating 'Lys-11'-linked polyubiquitin chains on APC/C substrates, leading to the degradation of APC/C substrates by the proteasome and promoting mitotic exit.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images