

BTN3A2 rabbit pAb

货号	YP-Ab-18004
同位型	IgG
应用	WB
种属	Human
靶点	BTN3A2
基因名称	BTN3A2 BT3.2 BTF3 BTF4
蛋白名称	Butyrophilin subfamily 3 member A2
免疫原	Synthesized peptide derived from human BTN3A2
特异性	This antibody detects endogenous levels of BTN3A2 at Human
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	37kD
功能	Plays a role in T-cell responses in the adaptive immune response. Inhibits the release of IFNG from activated T-cells.
细胞定位	Cell membrane ; Single-pass type I membrane protein .
组织表达	Detected in T-cells and natural killer cells.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images