



## MIDN rabbit pAb

货号	YP-Ab-18009
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	MIDN
基因名称	MIDN
蛋白名称	Midnolin (Midbrain nucleolar protein)
免疫原	Synthesized peptide derived from human MIDN
特异性	This antibody detects endogenous levels of MIDN at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	51kD
功能	Facilitates ubiquitin-independent proteasomal degradation of polycomb protein CBX4. Plays a role in inhibiting the activity of glucokinase GCK and both glucose-induced and basal insulin secretion.
细胞定位	Nucleus, nucleolus . Nucleus . Cytoplasm, cytosol . Detected in the nucleus and nucleolus with no expression in the cytoplasm (By similarity). However, a later study finds expression in the nucleus and cytoplasm with no expression in the nucleolus (PubMed:24187134). .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images