

ARMC10 rabbit pAb

货号	YP-Ab-18016
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	ARMC10
基因名称	ARMC10 SVH PSEC0198
蛋白名称	Armadillo repeat-containing protein 10 (Splicing variant involved in hepatocarcinogenesis protein)
免疫原	Synthesized peptide derived from human ARMC10
特异性	This antibody detects endogenous levels of ARMC10 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	38kD
功能	May play a role in cell survival and cell growth. May suppress the transcriptional activity of p53/TP53.
细胞定位	Endoplasmic reticulum membrane ; Single-pass membrane protein .
组织表达	Expressed in all tissues tested with higher expression in placenta, liver, kidney, heart and brain.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images