



SPANXB1 rabbit pAb

货号	YP-Ab-18021
同位型	IgG
应用	WB
种属	Human
靶点	SPANXB1
基因名称	SPANXB1 SPANXB; SPANXB2; SPANXF1
蛋白名称	Sperm protein associated with the nucleus on the X chromosome B/F (Cancer/testis antigen 11.2) (CT11.2) (Nuclear-associated protein SPAN-Xb) (SPANX-B) (Nuclear-associated protein SPAN-Xf) (SPANX-F)
免疫原	Synthesized peptide derived from human SPANXB1
特异性	This antibody detects endogenous levels of SPANXB1 at Human
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	11kD
细胞定位	Cytoplasm . Nucleus . Associated with nuclear craters.
组织表达	Detected in testis and sperm.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images