

MORC1 rabbit pAb

货号	YP-Ab-18023
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	MORC1
基因名称	MORC1 MORC
蛋白名称	MORC family CW-type zinc finger protein 1 (Cancer/testis antigen 33) (CT33)
免疫原	Synthesized peptide derived from human MORC1
特异性	This antibody detects endogenous levels of MORC1 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	108kD
细胞定位	Nucleus .
功能	Required for spermatogenesis (By similarity). Essential for de novo DNA methylation and silencing of transposable elements in the male embryonic germ cells (By similarity).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images