



SPATA19 rabbit pAb

货号	YP-Ab-18026
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	SPATA19
基因名称	SPATA19 SPERGEN1
蛋白名称	Spermatogenesis-associated protein 19, mitochondrial (Spermatogenic cell-specific gene 1 protein) (Spergen-1)
免疫原	Synthesized peptide derived from human SPATA19
特异性	This antibody detects endogenous levels of SPATA19 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	18kD
细胞定位	Mitochondrion outer membrane . Mitochondrion . Cell projection, cilium, flagellum . Localizes to the midpiece of the sperm flagellum. .
功能	Essential for sperm motility and male fertility (By similarity). Plays an important role in sperm motility by regulating the organization and function of the mitochondria and is also required for correct sperm midpiece assembly (By similarity).
组织表达	Expressed specifically in testis.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images