

NXF2 rabbit pAb

货号	YP-Ab-18028
同位型	IgG
应用	WB
种属	Human
靶点	NXF2; NXF2B
基因名称	NXF2 TAPL2; NXF2B
蛋白名称	Nuclear RNA export factor 2 (Cancer/testis antigen 39) (CT39) (TAP-like protein 2) (TAPL-2)
免疫原	Synthesized peptide derived from human NXF2
特异性	This antibody detects endogenous levels of NXF2 at Human
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	69kD
功能	Involved in the export of mRNA from the nucleus to the cytoplasm.
细胞定位	Nucleus, nucleoplasm. Cytoplasm. Localized in the nucleoplasm and at the nuclear envelope. Shuttles between the nucleus and the cytoplasm.
组织表达	Expressed almost exclusively in testis. Also expressed in several cancers.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images