

ODF1 rabbit pAb

货号	YP-Ab-18031
同位型	IgG
应用	WB
种属	Human,Rat
靶点	ODF1
基因名称	ODF1 ODFP
蛋白名称	Outer dense fiber protein 1
免疫原	Synthesized peptide derived from human ODF1
特异性	This antibody detects endogenous levels of ODF1 at Human,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	28kD
功能	Component of the outer dense fibers (ODF) of spermatozoa. ODF are filamentous structures located on the outside of the axoneme in the midpiece and principal piece of the mammalian sperm tail and may help to maintain the passive elastic structures and elastic recoil of the sperm tail.
组织表达	Testis.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images