



GSTA2 rabbit pAb

货号	YP-Ab-18032
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	GSTA2
基因名称	GSTA2 GST2
蛋白名称	Glutathione S-transferase A2 (EC 2.5.1.18) (GST HA subunit 2) (GST class-alpha member 2) (GST-gamma) (GSTA2-2) (GTH2)
免疫原	Synthesized peptide derived from human GSTA2
特异性	This antibody detects endogenous levels of GSTA2 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	24kD
功能	Catalyzes the conjugation of glutathione to a large variety of electrophilic compounds.
细胞定位	Cytoplasm.
组织表达	Liver.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images