



ACSS3 rabbit pAb

货号	YP-Ab-18039
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	ACSS3
基因名称	ACSS3
蛋白名称	Acyl-CoA synthetase short-chain family member 3, mitochondrial (EC 6.2.1.1) (Acyl-CoA synthetase short-chain family member 3)
免疫原	Synthesized peptide derived from human ACSS3
特异性	This antibody detects endogenous levels of ACSS3 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
分子量	75kD
功能	Catalyzes the synthesis of acetyl-CoA from short-chain fatty acids . Propionate is the preferred substrate . Can utilize acetate and butyrate with a much lower affinity (By similarity).
细胞定位	Mitochondrion matrix .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images