



cadherin pan (ABT172R) rabbit mAb

货号	YP-Ab-18047
同位型	IgG
应用	IHC;WB; ELISA
种属	Human; Mouse; (predicted: Rat)
靶点	cadherin pan
简介	>>Rap1 signaling pathway;>>Apelin signaling pathway;>>Hippo signaling pathway;>>Cell adhesion molecules;>>Adherens junction;>>Bacterial invasion of epithelial cells;>>Pathways in cancer;>>Endometrial cancer;>>Thyroid cancer;>>Melanoma;>>Bladder cancer;>>Gastric cancer
基因名称	CDH1/CDH2/CDH4
蛋白名称	cadherin
免疫原	Synthesized peptide derived from human cadherin-1 AA range:800-882
特异性	This antibody detects endogenous levels of cadherin pan
组成	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
来源	Monoclonal, Rabbit IgG1, Kappa
稀释	IHC 1:100-500, WB 1:500-1000, ELISA 1:5000-20000
纯化工艺	Recombinant Expression and Affinity purified
其他名称	Cadherin-1;CAM 120/80;Epithelial cadherin;E-cadherin;Uvomorulin;CD324;Cadherin-2;CDw325;Neural cadherin;N-cadherin1 Publication;CD325;Cadherin-3;Placental cadherin;P-cadherin;Cadherin-4;Retinal cadherin;R-CAD;R-cadherin
细胞定位	Membranous
组织表达	Non-neural epithelial tissues.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

