



OXSM rabbit pAb

货号	YP-Ab-18048
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	OXSM
基因名称	OXSM
蛋白名称	3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial (EC 2.3.1.41) (Beta-ketoacyl-ACP synthase)
免疫原	Synthesized peptide derived from human OXSM
特异性	This antibody detects endogenous levels of OXSM at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	50kD
功能	May play a role in the biosynthesis of lipoic acid as well as longer chain fatty acids required for optimal mitochondrial function.
细胞定位	Mitochondrion .
组织表达	Widely expressed. Higher expression in heart, skeletal muscle, liver and kidney which contain high levels of active mitochondria.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images