



## UXT rabbit pAb

货号	YP-Ab-18051
同位型	IgG
应用	WB
种属	Human,Rat
靶点	UXT
基因名称	UXT HSPC024
蛋白名称	Protein UXT (Androgen receptor trapped clone 27 protein) (ART-27) (Ubiquitously expressed transcript protein)
免疫原	Synthesized peptide derived from human UXT
特异性	This antibody detects endogenous levels of UXT at Human,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	17kD
功能	Involved in gene transcription regulation . Acts in concert with the corepressor URI1 to regulate androgen receptor AR-mediated transcription . Together with URI1, associates with chromatin to the NKX3-1 promoter region . Negatively regulates the transcriptional activity of the estrogen receptor ESR1 by inducing its translocation into the cytoplasm . May act as nuclear chaperone that facilitates the formation of the NF-kappa-B enhanceosome and thus positively regulates NF-kappa-B transcription activity . Potential component of mitochondrial-associated LRPPRC, a multidomain organizer that potentially integrates mitochondria and the microtubular cytoskeleton with chromosome remodeling . Increasing concentrations of UXT contributes to progressive aggregation of mitochondria and cell death potentially through its association with LRPPRC . Suppresses cell transformation and it might mediate
细胞定位	[Isoform 1]: Cytoplasm .; Cytoplasm . Nucleus . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole . Predominantly localizes to the nucleus (PubMed:16221885). Localizes to spindle pole during mitosis (PubMed:16221885). .
组织表达	Ubiquitous (PubMed:10087202, PubMed:11854421, PubMed:17635584, PubMed:11827465). Expressed in prostate epithelial cells (PubMed:21730289). Expressed in mammary epithelial cells (PubMed:28106301). Highest levels in the heart, skeletal muscle, pancreas, kidney, liver, adrenal gland, peripheral blood leukocytes, lymph node, prostate, and thyroid and the lowest levels in bladder and uterus (PubMed:11854421, PubMed:17635584, PubMed:11827465). Overexpressed in a number of tumor tissues (PubMed:11854421, PubMed:16221885, PubMed:28106301).
浓度	1 mg/ml

储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images