

## FAM50A rabbit pAb

货号	YP-Ab-18055
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	FAM50A
基因名称	FAM50A DXS9928E HXC26 XAP5
蛋白名称	Protein FAM50A (Protein HXC-26) (Protein XAP-5)
免疫原	Synthesized peptide derived from human FAM50A
特异性	This antibody detects endogenous levels of FAM50A at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	37kD
功能	Probably involved in the regulation of pre-mRNA splicing.
细胞定位	Nucleus .
组织表达	Widely expressed in fetal and adult tissues. Mostly abundant in fetal brain, liver and kidney; in the adult, high levels were also observed in heart, skeletal muscle, spleen, thymus, prostate and small intestine. Expressed in fetal cerebellum and hypothalamus. Low expression is observed in fetal temporal lobe (PubMed:32703943).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images