

SPATC1 rabbit pAb

货号	YP-Ab-18057
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	SPATC1
基因名称	SPATC1 SPATA15 SPRN
蛋白名称	Speriolin (Spermatogenesis and centriole-associated protein 1) (Spermatogenesis-associated protein 15) (Spermatogenic cell-specific Cdc20-binding protein)
免疫原	Synthesized peptide derived from human SPATC1
特异性	This antibody detects endogenous levels of SPATC1 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	65kD
细胞定位	Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Colocalizes with the centrosomal pericentrin protein PCNT1 (By similarity). Located in the connecting piece of sperm. .
组织表达	Detected only in testis.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images