



VMA21 rabbit pAb

货号	YP-Ab-18061
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	VMA21
基因名称	VMA21 MEAX XMEA
蛋白名称	Vacuolar ATPase assembly integral membrane protein VMA21 (Myopathy with excessive autophagy protein)
免疫原	Synthesized peptide derived from human VMA21
特异性	This antibody detects endogenous levels of VMA21 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	11kD
功能	Required for the assembly of the V0 complex of the vacuolar ATPase (V-ATPase) in the endoplasmic reticulum.
细胞定位	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Endoplasmic reticulum-Golgi intermediate compartment membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, COPII-coated vesicle membrane ; Multi-pass membrane protein .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images