

## GPLOW rabbit pAb

货号	YP-Ab-18063
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	GPLOW
基因名称	GPLOW GPATC5 GPATCH5 T54
蛋白名称	G patch domain and LOW motifs-containing protein (G patch domain-containing protein 5) (Protein MOS2 homolog) (Protein T54)
免疫原	Synthesized peptide derived from human GPLOW
特异性	This antibody detects endogenous levels of GPLOW at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	52kD
功能	RNA-binding protein involved in pre-mRNA splicing.
细胞定位	Nucleus .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images