



## RHOC rabbit pAb

|        |   |
|--------|---|
| 货号     | YP-Ab-18065   |
| 同位型    | IgG   |
| 应用     | WB  |
| 种属     | Human;Mouse   |
| 靶点     | RHOC  |
| 基因名称   | RHOC ARH9 ARHC  |
| 蛋白名称   | Rho-related GTP-binding protein RhoC (Rho cDNA clone 9) (h9)  |
| 免疫原    | Synthesized peptide derived from human RHOC   |
| 特异性    | This antibody detects endogenous levels of RHOC at Human, Mouse   |
| 组成     | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| 来源     | Polyclonal, Rabbit,IgG  |
| 稀释     | WB 1:500-2000   |
| 纯化工艺   | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.   |
| 分子量    | 21kD  |
| 功能     | Regulates a signal transduction pathway linking plasma membrane receptors to the assembly of focal adhesions and actin stress fibers. Serves as a microtubule-dependent signal that is required for the myosin contractile ring formation during cell cycle cytokinesis. Regulates apical junction formation in bronchial epithelial cells. |
| 细胞定位   | Cell membrane ; Lipid-anchor ; Cytoplasmic side . Cleavage furrow . Translocates to the equatorial region before furrow formation in a ECT2-dependent manner.   |
| 浓度     | 1 mg/ml   |
| 储存     | -15°C to -25°C/1 year(Do not lower than -25°C)  |
| 有关注意事项 | Avoid repeated freezing and thawing!  |
| 使用建议   | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.   |

### Products Images