



C10orf46 rabbit pAb

货号	YP-Ab-18075
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	CACUL1
基因名称	CACUL1 C10orf46 CAC1
蛋白名称	CDK2-associated and cullin domain-containing protein 1 (Cdk-associated cullin1)
免疫原	Synthesized peptide derived from human C10orf46
特异性	This antibody detects endogenous levels of C10orf46 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	41kD
功能	Cell cycle associated protein capable of promoting cell proliferation through the activation of CDK2 at the G1/S phase transition.
组织表达	Ubiquitously expressed with highest expression in the mammary gland and large intestine. Highly expressed in cancer tissues and cancer cell lines. During cell cycle progression expression is high in G1/S, low in the middle of S phase, and increases again in S/G2 phase.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images