

# MTP18 rabbit pAb

货号	YP-Ab-18078
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	MTFP1
基因名称	MTFP1 MTP18 HSPC242 My022
蛋白名称	Mitochondrial fission process protein 1 (Mitochondrial 18 kDa protein) (MTP18)
免疫原	Synthesized peptide derived from human MTP18
特异性	This antibody detects endogenous levels of MTP18 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	18kD
功能	Involved in the mitochondrial division probably by regulating membrane fission. Loss-of-function induces the release of cytochrome c, which activates the caspase cascade and leads to apoptosis.
细胞定位	Mitochondrion inner membrane ; Multi-pass membrane protein .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images