

# SLC25A10 rabbit pAb

货号	YP-Ab-18081
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	SLC25A10
基因名称	SLC25A10 DIC
蛋白名称	Mitochondrial dicarboxylate carrier (Solute carrier family 25 member 10)
免疫原	Synthesized peptide derived from human SLC25A10
特异性	This antibody detects endogenous levels of SLC25A10 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	32kD
功能	Involved in translocation of malonate, malate and succinate in exchange for phosphate, sulfate, sulfite or thiosulfate across mitochondrial inner membrane.
细胞定位	Mitochondrion inner membrane; Multi-pass membrane protein.
组织表达	Present in high amounts in liver and kidney, and at lower levels in all the other tissues analyzed.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images