



## ETF A rabbit pAb

货号	YP-Ab-18083
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	ETF A
基因名称	ETF A
蛋白名称	Electron transfer flavoprotein subunit alpha, mitochondrial (Alpha-ETF)
免疫原	Synthesized peptide derived from human ETF A
特异性	This antibody detects endogenous levels of ETF A at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	37kD
功能	Heterodimeric electron transfer flavoprotein that accepts electrons from several mitochondrial dehydrogenases, including acyl-CoA dehydrogenases, glutaryl-CoA and sarcosine dehydrogenase . It transfers the electrons to the main mitochondrial respiratory chain via ETF-ubiquinone oxidoreductase (ETF dehydrogenase) . Required for normal mitochondrial fatty acid oxidation and normal amino acid metabolism .
细胞定位	Mitochondrion matrix .
浓度	1 mg/ml
储存	-15℃ to -25℃/1 year(Do not lower than -25℃)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images