



SMTN rabbit pAb

货号	YP-Ab-18091
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	SMTN
基因名称	SMTN SMSMO
蛋白名称	Smoothelin
免疫原	Synthesized peptide derived from human SMTN
特异性	This antibody detects endogenous levels of SMTN at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	101kD
功能	Structural protein of the cytoskeleton.
细胞定位	Cytoplasm, cytoskeleton. Exhibits a filamentous organization.
组织表达	Smooth muscle; contractile or vascular (for the long form).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images