



## PYGL rabbit pAb

货号	YP-Ab-18098
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	PYGL
基因名称	PYGL
蛋白名称	Glycogen phosphorylase, liver form (EC 2.4.1.1)
免疫原	Synthesized peptide derived from human PYGL
特异性	This antibody detects endogenous levels of PYGL at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	93kD
功能	Allosteric enzyme that catalyzes the rate-limiting step in glycogen catabolism, the phosphorolytic cleavage of glycogen to produce glucose-1-phosphate, and plays a central role in maintaining cellular and organismal glucose homeostasis.
细胞定位	Cytoplasm, cytosol .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images