

SIDT1 rabbit pAb

货号	YP-Ab-18104
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	SIDT1
基因名称	SIDT1
蛋白名称	SID1 transmembrane family member 1
免疫原	Synthesized peptide derived from human SIDT1
特异性	This antibody detects endogenous levels of SIDT1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	91kD
功能	In vitro binds long double-stranded RNA (dsRNA) (500 and 700 base pairs), but not dsRNA shorter than 300 bp. Not involved in RNA autophagy, a process in which RNA is directly imported into lysosomes in an ATP-dependent manner, and degraded.
细胞定位	Membrane ; Multi-pass membrane protein .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images