



## SLC38A5 rabbit pAb

货号	YP-Ab-18108
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	SLC38A5
基因名称	SLC38A5 JM24 SN2 SNAT5 PP7194
蛋白名称	Sodium-coupled neutral amino acid transporter 5 (Solute carrier family 38 member 5) (System N transporter 2)
免疫原	Synthesized peptide derived from human SLC38A5
特异性	This antibody detects endogenous levels of SLC38A5 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	52kD
功能	Functions as a sodium-dependent amino acid transporter which countertransport protons. Mediates the saturable, pH-sensitive, and electrogenic cotransport of several neutral amino acids including glycine, asparagine, alanine, serine, glutamine and histidine with sodium.
细胞定位	Cell membrane ; Multi-pass membrane protein .
组织表达	Predominantly expressed in stomach, brain, liver, lung and intestinal tract.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

