

FIBIN rabbit pAb

货号	YP-Ab-18109
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	FIBIN
基因名称	FIBIN PSEC0235
蛋白名称	Fin bud initiation factor homolog
免疫原	Synthesized peptide derived from human FIBIN
特异性	This antibody detects endogenous levels of FIBIN at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	23kD
细胞定位	Secreted . Golgi apparatus . Endoplasmic reticulum
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images