



TDO2 rabbit pAb

货号	YP-Ab-18119
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	TDO2
基因名称	TDO2 TDO
蛋白名称	Tryptophan 2,3-dioxygenase (TDO) (EC 1.13.11.11) (Tryptamin 2,3-dioxygenase) (Tryptophan oxygenase) (TO) (TRPO) (Tryptophan pyrrolase) (Tryptophanase)
免疫原	Synthesized peptide derived from human TDO2
特异性	This antibody detects endogenous levels of TDO2 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	45kD
功能	Heme-dependent dioxygenase that catalyzes the oxidative cleavage of the L-tryptophan (L-Trp) pyrrole ring and converts L-tryptophan to N-formyl-L-kynurenine. Catalyzes the oxidative cleavage of the indole moiety.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images